



Même®



Même® Microthane®

kulatá base
kulatý profil
nízká projekce (LP)


Katalog.číslo	A/B [mm]	C [mm]	D [mm]	Objem [ml]
30724-055*1	73	25	46	55
30724-065*1	78	26	49	65
30724-085*1	83	26	51	85
30724-095*1	88	27	54	95
30724-110*1	93	28	57	110
30724-125*1	98	29	60	125
30724-150*1	103	30	62	150
30724-165*1	108	31	65	165
30724-185*	113	31	68	185
30724-215	118	32	71	215
30724-235*	123	33	73	235
30724-270	128	34	76	270
30724-295*1	133	35	79	295
30724-330*1	138	36	82	330
30724-365*1	143	37	85	365
30724-400*1	148	37	87	400
30724-440*1	153	38	90	440
30724-485*1	158	39	93	485



Même® Microthane®

kulatá base
kulatý profil
střední projekce (MP)


Katalog.číslo	A/B [mm]	C [mm]	D [mm]	Objem [ml]
30725-070*1	73	30	51	70
30725-090*1	78	32	54	90
30725-105*1	83	33	57	105
30725-120*1	88	34	60	120
30725-145*1	93	35	63	145
30725-165*1	98	37	66	165
30725-195*	103	38	69	195
30725-220	108	39	72	220
30725-255*	113	40	75	255
30725-285	118	42	78	285
30725-320*	123	43	81	320
30725-360	128	45	84	360
30725-400*	133	46	87	400
30725-445	138	47	90	445
30725-495*	143	49	93	495
30725-545*1	148	50	96	545
30725-600*1	153	51	99	600
30725-660*1	158	53	102	660



Même® Microthane®

kulatá base
kulatý profil
vysoká projekce (HP)

Katalog.číslo	A/B [mm]	C [mm]	D [mm]	Objem [ml]
30726-090*1	73	36	55	90
30726-110*1	78	38	58	110
30726-125*1	83	39	62	125
30726-155*1	88	41	65	155
30726-180*	93	42	68	180
30726-210	98	44	71	210
30726-235*	103	46	75	235
30726-275	108	47	78	275
30726-315*	113	49	81	315
30726-360	118	51	84	360
30726-400*	123	53	87	400
30726-450	128	55	91	450
30726-505*	133	57	94	505
30726-560	138	59	97	560
30726-620*1	143	61	100	620
30726-690*1	148	63	104	690
30726-760*1	153	64	107	760
30726-835*1	158	66	110	835




Même® Microthane®

kulatá base
kulatý profil
extra vysoká projekce (XP)

Katalog.číslo	A/B [mm]	C [mm]	D [mm]	Objem [ml]
30727-105*1	73	42	59	105
30727-125*1	78	45	63	125
30727-155*1	83	47	66	155
30727-180*1	88	49	70	180
30727-215*1	93	51	73	215
30727-250*1	98	53	77	250
30727-285*1	103	55	81	285
30727-330	108	57	84	330
30727-380*	113	59	88	380
30727-430	118	61	91	430
30727-485*	123	63	95	485
30727-545	128	65	98	545
30727-610*	133	67	102	610
30727-675	138	69	105	675
30727-785*1	143	71	109	755
30727-830*1	148	73	112	830
30727-920*1	153	75	116	920
30727-1010*1	158	77	119	1010


Replicon®



Replicon® Microthane®

kulatá base
anatomický profil
nízká projekce (LP)


Katalog.číslo	A/B [mm]	C [mm]	D [mm]	Objem [ml]
30734-050*1	73	25	34	50
30734-060*1	78	26	36	60
30734-070*1	83	26	38	70
30734-085*1	88	27	40	85
30734-100*1	93	28	42	100
30734-110*1	98	29	44	110
30734-125*1	103	30	46	125
30734-150*1	108	31	48	150
30734-165*	113	31	50	165
30734-185	118	32	52	185
30734-210*	123	33	54	210
30734-235	128	34	56	235
30734-265*	133	35	58	265
30734-290	138	36	59	290
30734-320*	143	37	61	320
30734-350*1	148	37	63	350
30734-390*1	153	38	65	390
30734-425*1	158	39	67	425



Replicon® Microthane®

kulatá base
anatomický profil
střední projekce (MP)


Katalog.číslo	A/B [mm]	C [mm]	D [mm]	Objem [ml]
30735-065*1	73	30	41	65
30735-075*1	78	32	43	75
30735-095*1	83	33	45	95
30735-110*1	88	34	48	110
30735-125*1	93	35	50	125
30735-150	98	37	52	150
30735-170*	103	38	54	170
30735-195	108	39	56	195
30735-220*	113	40	59	220
30735-255	118	42	61	255
30735-280*	123	43	63	280
30735-320	128	45	65	320
30735-350*	133	46	67	350
30735-390	138	47	70	390
30735-435*	143	49	72	435
30735-480	148	50	74	480
30735-530*	153	51	76	530
30735-580*1	158	53	78	580



Replicon® Microthane®

kulatá base
anatomický profil
vysoká projekce (HP)

Katalog.číslo	A/B [mm]	C [mm]	D [mm]	Objem [ml]
30736-075*1	73	36	46	75
30736-095*1	78	38	48	95
30736-115*1	83	39	51	115
30736-130*1	88	41	53	130
30736-160*1	93	42	56	160
30736-180	98	44	58	180
30736-210*	103	46	61	210
30736-240	108	47	63	240
30736-275*	113	49	66	275
30736-315	118	51	68	315
30736-350*	123	53	71	350
30736-395	128	55	73	395
30736-445*	133	57	76	445
30736-495	138	59	78	495
30736-550*	143	61	81	550
30736-605	148	63	83	605
30736-670*1	153	64	86	670
30736-735*1	158	66	88	735



Replicon® Microthane®

kulatá base
anatomický profil
extra vysoká projekce (XP)

Katalog.číslo	A/B [mm]	C [mm]	D [mm]	Objem [ml]
30737-095*1	73	42	50	95
30737-110*1	78	45	53	110
30737-130*1	83	47	56	130
30737-160*1	88	49	59	160
30737-185*1	93	51	62	185
30737-220	98	53	65	220
30737-255*	103	55	67	255
30737-290	108	57	70	290
30737-330*	113	59	73	330
30737-380	118	61	76	380
30737-425*	123	63	79	425
30737-480	128	65	82	480
30737-535*	133	67	84	535
30737-600	138	69	87	600
30737-665*	143	71	90	665
30737-730*1	148	73	93	730
30737-810*1	153	75	96	810
30737-890*1	158	77	99	890

* Dodávka po objednání = 6 týdnů

* Sizer skladem; Katalog.číslo e.g. e.g. 20725-195H

